

STORAGE AND ACCESSORIES

STORAGE PRODUCTS AND AUTOMATION ACCESSORIES



PAA PROVIDES AN EXTENSIVE RANGE OF AUTOMATION ACCESSORIES REQUIRED TO COMPLETE A FULLY AUTOMATED SOLUTION.

THESE ARE DESIGNED TO BE INTEGRATED BOTH WITH OUR RANGE OF LABORATORY ROBOTS AND OTHER UNIVERSALLY USED ROBOTIC SYSTEMS.

Hotel rack storage

Hotel Racks are designed to store microplates for quick, at-will retrieval to load/unload media. Each PAA Hotel Rack is compatible with various PAA Docking Station and Carousel products.

Hotel rack universal

PAA's Hotel Rack Universal (HRU) provides:

- / **Faster processing**
It stores microplates in both landscape and portrait orientation, eliminating the need for re-grip nests
- / **Flexible solution**
Each shelf can be removed to create additional spacing for tip racks or deep well plates
- / **Adaptable capacity**
12, 17 or 20 universal shelves
- / **Compliance**
RoHS and CE marked.

Sequential stack storage

Sequential Stacks offer rapid retrieval to load/unload media and quickly return it back to the assay.

- / **Flexible capacity**
PAA offers an extensive range of Sequential Stacks with capacity for 30, 45 and 50 SBS microplates
- / **Choice of orientation**
Available in landscape or portrait orientation
- / **Compliance**
RoHS and CE marked.



AT THE HEART
OF LABORATORY
AUTOMATION

Docking Stations

PAA's Docking Stations are intelligently engineered to repeatedly capture and release any configuration of PAA Hotel Racks or Sequential Stacks.



Docking Station Universal

/ Simple installation

PAA's Docking Station Universal (DSU) offers adjustable levelling feet and surface mount locations

/ Ease of use

A quick release button at rear of the docking station releases the storage device for easy retrieval

/ Adaptable capacity

The attachment bracket allows for easy modular connection with other DSUs for increased capacity

/ Total Flexibility

DSUs can be assembled in multiple radius configurations tailored to your requirements

/ Compliance

RoHS and CE marked.

Automation accessories



Barcode Reader




Nests



Teach Plates


Optical mount plates:



GX0P



KX0P



CS0P

PAA

AT THE HEART
OF LABORATORY
AUTOMATION

United Kingdom
Peak Analysis and Automation Ltd
T: +44 (0)1252 373000

North America
Peak Analysis & Automation Inc
T: +1-719-598-3555

info@paa-automation.com paa-automation.com